


# Fixture Size BD Cuff IMPLANT

## NC Ø3.7

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c3.2)Ø3.7	NC	5	2	7	Ø4.0	Ø4.4	BDC370502C	BDC370502DC
			4	9			BDC370504C	BDC370504DC
			6	11			BDC370506C	BDC370506DC
			8	13			BDC370508C	BDC370508DC
			10	15			BDC370510C	BDC370510DC
		7	2	9			BDC370702C	BDC370702DC
			4	11			BDC370704C	BDC370704DC
			6	13			BDC370706C	BDC370706DC
			8	15			BDC370708C	BDC370708DC
		9	2	11			BDC370902C	BDC370902DC
			4	13			BDC370904C	BDC370904DC
			6	15			BDC370906C	BDC370906DC



## RC Ø4.1

- Cover Screw 포함

Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c3.5)Ø4.1	RC	5	2	7	Ø4.4	Ø4.8	BDC410502C	BDC440502DC
			4	9			BDC410504C	BDC440504DC
			6	11			BDC410506C	BDC440506DC
			8	13			BDC410508C	BDC410508DC
			10	15			BDC410510C	BDC410510DC
		7	2	9			BDC410702C	BDC440702DC
			4	11			BDC410704C	BDC440704DC
			6	13			BDC410706C	BDC440706DC
			8	15			BDC410708C	BDC410708DC
		9	2	11			BDC410902C	BDC440902DC
			4	13			BDC410904C	BDC440904DC
			6	15			BDC410906C	BDC440906DC




# Fixture Size BD Cuff IMPLANT

## RC Ø4.4

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c3.8)Ø4.4	RC	5	2	7	Ø4.7	Ø5.1	BDC440502C	BDC440502DC
			4	9			BDC440504C	BDC440504DC
			6	11			BDC440506C	BDC440506DC
			8	13			BDC410508C	BDC440508DC
			10	15			BDC410510C	BDC440510DC
		7	2	9			BDC440702C	BDC440702DC
			4	11			BDC440704C	BDC440704DC
			6	13			BDC440706C	BDC440706DC
			8	15			BDC410708C	BDC440708DC
			2	11			BDC440902C	BDC440902DC
		9	4	13			BDC440904C	BDC440904DC
			6	15			BDC440906C	BDC440906DC



## RC Ø4.8

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c4.2)Ø4.8	RC	5	2	7	Ø4.7	Ø5.1	BDC480502C	BDC480502DC
			4	9			BDC480504C	BDC480504DC
			6	11			BDC480506C	BDC480506DC
		7	2	9			BDC480702C	BDC480702DC
			4	11			BDC480704C	BDC480704DC
			6	13			BDC480706C	BDC480706DC
		9	2	11			BDC480902C	BDC480902DC
			4	13			BDC480904C	BDC480904DC
			6	15			BDC480906C	BDC480906DC



## RC Ø5.3

- Cover Screw 포함

Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c4.5)Ø5.3	RC	5	2	7	Ø5.6	Ø6.0	BDC530502C	BDC530502DC
			4	9			BDC530504C	BDC530504DC
		7	2	9			BDC530702C	BDC530702DC
			4	11			BDC530704C	BDC530704DC




# Fixture Size BD Cuff IMPLANT

## RC Ø5.8

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c4.8)Ø5.8	RC	5	2	7	-	Ø6.4	-	BDC580502DC
			4	9			-	BDC580504DC
		7	2	9			-	BDC580702DC
			4	11			-	BDC580704DC



## RC Ø6.3

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c4.8)Ø6.3	RC	5	2	7	-	Ø6.9	-	BDC630502DC
			4	9			-	BDC630504DC
		7	2	9			-	BDC630702DC
			4	11			-	BDC630704DC



## RC Ø6.8

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c4.8)Ø6.8	RC	5	2	7	-	Ø7.4	-	BDC680502DC
			4	9			-	BDC680504DC
		7	2	9			-	BDC680702DC
			4	11			-	BDC680704DC



## RC Ø7.3

- Cover Screw 포함

Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Total Length (mm)	Widest thread diameter		Ref.C	
					Normal	Deep	Normal	Deep
(c4.8)Ø7.3	RC	5	2	7	-	Ø7.9	-	BDC730502DC
			4	9			-	BDC730504DC
		7	2	9			-	BDC730702DC
			4	11			-	BDC730704DC




# Fixture Size BD Cuff IMPLANT

## RC Ø5.3 Deep Thread

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Length (mm)	Ref.C
Ø5.3	RC	5	2	7	BDC530502DC
			4	9	BDC530504DC
		7	2	9	BDC530702DC
			4	11	BDC530704DC



## RC Ø5.8 Deep Thread

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Length (mm)	Ref.C
Ø5.8	RC	5	2	7	BDC580502DC
			4	9	BDC580504DC
		7	2	9	BDC580702DC
			4	11	BDC580704DC



## RC Ø6.3 Deep Thread

- Cover Screw 포함


Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Length (mm)	Ref.C
Ø6.3	RC	5	2	7	BDC630502DC
			4	9	BDC630504DC
		7	2	9	BDC630702DC
			4	11	BDC630704DC



## RC Ø6.8 Deep Thread

- Cover Screw 포함

Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Length (mm)	Ref.C
Ø6.8	RC	5	2	7	BDC680502DC
			4	9	BDC680504DC
		7	2	9	BDC680702DC
			4	11	BDC680704DC



## RC Ø7.3 Deep Thread

- Cover Screw 포함

Fixture Diameter	Connection	Thread Length (mm)	Cuff Height (mm)	Length (mm)	Ref.C
Ø7.3	RC	5	2	7	BDC730502DC
			4	9	BDC730504DC
		7	2	9	BDC730702DC
			4	11	BDC730704DC

